

NAME: _____
COMPANY/ORGANIZATION: _____
DEALER/REP: _____

PHONE: _____
E-MAIL: _____
SALES REP: _____ **QUOTE:** _____



PHSChiropractic
 by Pivotal Health Solutions

www.phschiropractic.com • 800-743-7738

ATT-300-W

Wood Roller Massage Table

STANDARD FEATURES:

- Load Capacity: 350 lbs.
- 72”L x 25”W x 25”H
- Hardwood construction
- 25” Wide top
- Clockwise and counter-clockwise rotation
- Travel indicator
- Automatic lowering device
- Adjustable massage pressure
- Vibration
- LED roller height indicator
- 30 Minute timer
- (1) Cervical bolster (2.5” x 12” x 14”) and (1) knee bolster (6” x 7” x 20”)
- 4” Roller height adjustment



ATT-300-W

Specially designed for pretreatment therapy and massage, the ATT-300-W Wood Roller Massage Table features three large rotating rollers traveling both clockwise and counterclockwise, up and down the spine, to knead and relax muscles in the neck, shoulders and lower back. For added patient comfort, the table features a 25” wide, thick padded lounge top, vibration, and adjustable massage pressure. A cervical and knee bolster are also included.

STANDARD TABLE:

- ATT-300-W** - Wood Roller Massage Table
 72”L x 25”W x 25”H

Qty _____

OPTIONS:

- ATT-RSC** - Reverse Side Controls

STANDARD UPHOLSTERY CHOICES:

- | | |
|--|-----------------------------------|
| <input type="checkbox"/> Imperial Blue | <input type="checkbox"/> Lark |
| <input type="checkbox"/> Blue Ridge | <input type="checkbox"/> Burgundy |
| <input type="checkbox"/> Putty | <input type="checkbox"/> Dove |
| <input type="checkbox"/> Forest | <input type="checkbox"/> Charcoal |
| <input type="checkbox"/> Black | <input type="checkbox"/> Mocha |

STANDARD WOOD CHOICES:

- | | |
|------------------------------|--------------------------------|
| <input type="checkbox"/> Oak | <input type="checkbox"/> Maple |
|------------------------------|--------------------------------|

WOOD STAIN OPTION:

- ATT-300-W-STAIN**

STAIN CHOICES:

- | | |
|--------------------------------------|--|
| <input type="checkbox"/> Warm Copper | <input type="checkbox"/> Rustic Fall |
| <input type="checkbox"/> Deep Garnet | <input type="checkbox"/> Warm Caramel |
| <input type="checkbox"/> Tungsten | <input type="checkbox"/> Custom Wood Stain Color |



ATT-300-W with Stain