## Chattanooga®

# INTELECT® RPW 2

THE NEXT GENERATION IN SHOCKWAVE THERAPY

## KEY FEATURES

- OPTIMAL ERGONOMIC DESIGN
- COMFORT MODE
- EASY FIT IN EVERY TREATMENT ROOM
- SUGGESTED PROTOCOLS LIBRARY
- LARGE GRAPHIC COLOR TOUCH SCREEN
- NEW FALCON® HANDPIECE
- NEW V-ACTOR® HF HANDPIECE (option)
- BLUETOOTH CONNECTIVITY WITH WINDOWS 10 APP
- EXTENDED PRESSURE RANGE from 0.3 to 5 bar
- LARGE CHOICE OF TRANSMITTERS (option)
- LARGE SELECTION OF LANGUAGES ONBOARD

## INDICATIONS

- ACHILLES TEHDINOPATHY (1)
- DISORDERS OF TENDONS INSERTIONS (2)
- MYOFASCIAL TRIGGER POINTS (3)
- PAIN & INFLAMMATION IN ORTHOPAEDIC CONDITIONS (4)
- PLANTAR FASCIITIS (5)
- PULSE VIBRATION MASSAGE (6)

### **TECHNICAL INFORMATION**

2176 KIT: 120V AC 50-60Hz 350VA

DIMENSIONS WxHxD: 16.15in x 16.54in x 43.31in WEIGHT: 66.14 Lbs COMPRESSED AIR OUTPUT: 0.3-5.0 bar POWER INCREMENT SETTINGS: 0.1 bar PULSE AMPLITUDE: Approx 2mm at 3 bar, idling (without coupling) PULSE WIDTH: Approx 6-12ms PULSE FREQUENCY: FALCON® handpiece 1-21Hz for RPW, V-ACTOR® handpiece 1-50Hz for vibration MODE: Single pulse or continuous USER INTERFACE: 7 inch, 800x480 pixels color graphic capacitive touch screen "CONFORMITY: This device complies with the applicable standards EN/IEC 60601-1, CAN/CSA-C22.2 No.601.1, UL Std. No 60601-1. CLASSIFICATION ACCORDING TO MDD: class IIa device 93/42 60601-1 LEGAL MANUFACTURER: DJO France SAS - CE 2797 chattanooga

# chattanooga® INTELECT® **RPW 2**

#### **ORDERING INFORMATION**

PART NUMBER	DESCRIPTION
2176KIT	Intelect® RPW 2 120V AC~50 - 60 Hz 350VA

### STANDARD ACCESSORIES

PART NUMBER	DESCRIPTION				
18080	FALCON $^{\circ}$ handpiece, including 2x projectile; 2x guiding tube; 1 cleaning brush, 1x wrench; 1x allen key				
1x Power cord with EU plug for 2174KIT, 1x Power cord with North American plug for 2173KIT, 1x Power cord with Japanese plug for 2175KIT					
13-28660*	Quick start guide booklets (available in multiple languages)	28178	Transmitter Ro40		
		29724	Transmitter D20-S black		
22652	Transmitters tray				
22651		82-0274	Intelect RPW 2 pouch		
22651	Gel bottle support cup	15-1140	USB Key		
4248	1x Transmission gel. 250 ml				

### **OPTIONAL ACCESSORIES**

PART NUMBER	DESCRIPTION			
81-0133	Intelect® RPW2 PRO Package-Optional Accessory Kit includes 19365 V-Actor HF handpiece with V25 & V40 Transmitters, the 29728 Di15 Gold Transmitter & 29539 Atlas (soft) RPW Transmitter.			
18080	FALCON handpiece, including 2x projectile; 2x guiding tube; 1 cleaning brush, 1x wrench; 1x allen key			
29104	Revision kit short including 2x projectiles; 2x guiding tubes short; 2x o-ring set; 1x o-ring guide			
29801	A6 transmitter black	82-0274	Intelect RPW 2 pouch	
29802	T10 transmitter black	22651	Gel bottle support cup	
29729	C15 transmitter black	22652	Transmitter tray	
29728	DI15 transmitter black	4248	24x Transmission gel, 250 ml	
29726	F15 transmitter black	14820	Power cord EU	
29724	D20-S transmitter black	14679	Power cord North America	
29725	D35-S transmitter black	14821	Power cord JAPAN	
28736	D20-T transmitter black	14814	Power cord Swiss	
28178	Ro40 transmitter red	14819	Power cord Australia	
18638	R15 transmitter black	14815	Power cord UK	
29539	Atlas RPW transmitter	14816	Power cord India	
19365	V-ACTOR® HF handpiece set Including V-ACTOR®, V25 and V40 transmitters	14817	Power cord Israel	
		14818	Power cord Denmark	
29742	V-ACTOR V10 transmitter			
28740	V-ACTOR V25 transmitter			
28741	V-ACTOR V40 transmitter			
31800	Transmitter box including 9 transmitters (A6, T10, C15, D115, F15, D20-S, D35-S, D20-T, Ro40)			
28945	PERI-ACTOR <sup>®</sup> fascia transmitter set including 4 transmitters			
28946	SPINE-ACTOR® transmitter set including 3 transmitters			

Falcon®, V-ACTOR®, SPINE-ACTOR®, and PERI-ACTOR® are registered trademarks of Storz Medical AC

1. Gerdesmeyer L. Current evidence of extracorporeal shock wave therapy in chronic Achilles tendinopathy NCBINCBI Logo Int J Surg. 2015 Dec;24(Pt B):154-9. doi: 10.1016/j ijsu. 2015.07.718. Epub 2015 Aug 29.\*

2. Gerdesmeuer L. Current evidence of extracorporeal shock wave therapy in chronic Achilles tendinopathy, NCBINCBI Logo Int J Surg, 2015

3. Nedelka T.Mechano-transduction effect of shockwaves in the treatment of lumbar facet joint pain.Neuro Endocrinol Lett. 2014

 Schmitz C. Efficacy and safety of extracorporeal shock wave therapy for orthopedic conditions: a systematic review on studies listed in the PEDro database. Br Med Bull. 2015;116:115-38. doi: 10.1093/bmb/ldv047. Epub 2015 Nov 18.

5. Ogden JA. Shockwave therapy for chronic proximal plantar fasciitis: a meta-analysis.Foot Ankle Int. 2002 Apr;23(4):301-8

6. S. Broadbent, JJ. Rousseau, RM. Thorp, SL. Choate, FS. Jackson, DS. Rowlands. Vibration therapy reduces plasma IL6 and muscle soreness after downhill running. 201 Sen:44(12):R8R-94

Individual results may vary. Neither DJO, LLC nor any of its subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any o type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.



T 800.336.6569 F 800.457.4221 1430 Decision Street I Vista, CA 92081-8553 I U.S.A. djoglobal.com/chattanooga

Copyright © 2021 by DJO, LLC MKT000-11172 Rev A