

# RPW Transmitters

A choice of standard transmitters for Radial Pressure Wave (RPW) indications. Each transmitter is unique by shape and material to deliver its own energy flux density and wave profile.



**TRIGGER POINTS/TENDINOPATHIES**

**R040 RED**



15 mm Energy beam transmitter with concave coupling surface, best for pain zones near the skin surface

Penetration depth 0-35mm

Part Number: 28178

**C15 BLACK**



Ceramic Energy 15mm transmitter for any type of tendonopathies

Penetration depth 0-35mm

Part Number: 29729

**D115 BLACK**



Golden Depth 15mm Deep Impact® transmitter for deep target areas, chronic disorders, local trigger points

Penetration depth 0-60mm

Part Number: 29728

**A6 BLACK**



Acupuncture like transmitter 6mm for local trigger points

Penetration depth 0-20mm

Part Number: 29801

**T10 BLACK**



Finger transmitter 10mm for deep target tendonopathies

Penetration depth 0-40mm

Part Number: 29802

**TRIGGER POINTS/TENDINOPATHIES**

**F15 BLACK**



White soft transmitter 15mm for superficial pain regions, muscles of mastication, cervical spine

Penetration depth 0-30mm

Part Number: 29726

**B15 BLACK**



White soft transmitter 15mm for superficial pain regions, muscles of mastication, cervical spine

Penetration depth 0-30mm

Part Number: 25377

**MYOFASCIAL TREATMENT**

**D20-S BLACK**




Standard Oscillator, 20mm transmitter or muscle and connective tissue

Penetration depth 0-50mm

Part Number: 29724

**D20-T BLACK**



Golden Oscillator 20mm transmitter for professional athletes, High energy application

Penetration depth 0-50mm

Part Number: 28736

**D35-S BLACK**



Large oscillator 35mm for large muscles and connective tissue

Penetration depth 0-60mm

Part Number: 29725

**VIBRATION THERAPY**

**V-ACTOR V10**



10mm for superficia trigger points

Part Number: 29742

**V-ACTOR V25**



25mm for muscle massage

Part Number: 28740

**V-ACTOR V40**



40mm for large muscles groups

Part Number: 28741

**TRIGGER POINTS/TENDINOPATHIES**

**ATLAS RPW**



Soft flexible tipped transmitter 15mm for highly sensitive regions and triggerpoints in the cervical spine

Penetration depth 0-15mm

Part Number: 29539

**SPINE TREATMENTS**

**SPINE-ACTORS®**



Double Tip shape to allows treatments of the paravertebral region.

3 sizes span from 20 to 35mm

Part Number: 28946

**FASCIA TREATMENTS**

**PERI-ACTORS®**



Knuckle



Scoop



Scraper



Sphere

Unique shape especially designed to dissolve clots in the Fascia structures

**PERI-ACTOR® Knuckle:** For Medium-sized and large fasciae thoracolumbar fasciae, plantar fasciae, trigger bands, shoulder and neck region.

**PERI-ACTOR® Scoop:** Treatment at or below the edges of bones: Front iliac crest, edge of the tibia

**PERI-ACTOR® Scraper:** Treatment of all fasciae trigger band

**PERI-ACTOR® Sphere:** Punctual treatment of fasciae: Plantar fasciae, Arch of foot, Hand, Trigger points, Cervical spine Fasciae

Part Number: 28945