

DOC DECOMPRESSION TABLE CHART

Is the DOC Decompression Table Right For Your Practice?

The unique design of the **DOC Decompression table** utilizes a state-of-the-art digital command center for specific vertebral targeting, separate lumbar and cervical traction programming and continuous readout and graphing of treatment protocols. In addition, the digital command center also controls table elevation at the push of a button.



Treats Prone & Supine
Treats in Flexion & Extension

The DOC Decompression Table positions so you can treat...



...in lateral flexion



...in rotation

Continually Monitors Feedback and Fully Customizable Protocols

The 8-inch, full color screen has a built-in Disc Angle Level Indicator that allows you to quickly target specific vertebral discs, making treatments accurate and effective.



Separate Lumbar and Cervical Motors

Instead of using the same traction headpiece to treat both the lumbar and cervical spine (which requires spinning the patient 180 degrees and puts the patients head where the previous patients' feet were), the DOC has separate lumbar and cervical motors that operate independently, so there is no need to change patient position when switching from lumbar to cervical. This saves time and increases efficiency.

No ropes or pulleys...the table moves to give more accurate pull for traction.



Belt Driven Decompression
Full Lumbar Harness

- ✓ Ships Ready to Treat
- ✓ Training & Marketing Available
- ✓ FDA Cleared
- ✓ CE Cleared for Canada
- ✓ TGA Cleared for Australia

Pivotal
PROFESSIONAL



Dr. Brian Self
DOC Decompression

A marketing program and technique training program is also available through Dr. Brian Self with Decompression Pros (not included).



800-743-7738
www.PHSChiropractic.com