

NAME: _____

PHONE: _____

COMPANY/ORGANIZATION: _____

E-MAIL: _____

DEALER/REP: _____

SALES REP: _____ QUOTE: _____



PHSChiropractic
by Pivotal Health Solutions

www.phschiropractic.com • 800-743-7738

TradeFlex Manual Flexion Table



STANDARD FEATURES:

- Stationary table height is 20"
- P/A Flexion 15°
- Rear lateral flexion 15°
- Electric spring tension
- 10" Ankle extension
- Breakaway abdominal
- 22" Wide comfort cushions
- Tilting headpiece with adjustable face cushions
- Black base
- Short T bar and SST strap
- Paper roller, cutter and guide bar

The TradeFlex Manual Flexion table is a workhorse, giving you and your patients a wide variety of options, from drops and headpieces to height, elevation and even long-axis distraction. You will appreciate the electric tension controls and removable flexion control t-bar that will allow you to move freely around the table. Your patients will benefit from the ability to perform figure-8 range of motion traction, a smooth range of 15 degrees of lateral flexion and a flexion depth of 15 degrees to accommodate almost any angle of drive.

The TradeFlex Manual Flexion table also offers versatility that will save you time and money, allowing for optimum correction of spinal subluxations or disc problems yet still functioning as an exam, adjusting or therapy table.

SPECIFICATIONS:

- 350 lbs. lifting capacity
- Overall dimensions: 73"L x 28"W x 20"H
- 83" Extended length (with ankle extension)
- Cushion width: 22"
- Rated voltage: 120V/230V
- Rated frequency: 60Hz

STANDARD TABLES:

- E9017** - TradeFlex Manual Flexion Table (120V)

Qty _____

- E9217** - TradeFlex Manual Flexion Table (230V)

Qty _____

HEIGHT & ELEVATION OPTIONS:

- E4802** - Extended stationary height beyond 26" (per inch) (Stationary table height is 20" - Adding height is no extra charge up to 26")

Qty _____ (inches)

- E9017-L** - Elevation - 21-29" with single foot pedal

- 95284** - Dual foot pedal for elevation



TradeFlex

Manual Flexion Table

HEADPIECE OPTIONS:

- E4401** - Tilt, Forward Drop, Adjustable Face Cushion
- E4408** - Tilt, Elevation, Adjustable Face Cushion
- E4409** - Tilt, Cervical Distraction, Adjustable Face Cushion
- E4410** - Tilt, Elevation, Forward Drop, Toggle Drop, Adjustable Face Cushion
- E4411** - Tilt, Elevation, Forward Drop, Toggle Drop, Cervical Distraction, Adjustable Face Cushion
- E4423** - Lateral Flexion (14") - Available on all headpiece options

LONG AXIS DISTRACTION OPTIONS:

- E4415** - Long axis distraction (120V)
- E4416** - Long axis distraction (230V)

DROP OPTIONS:

- E4412** - Thoracic Manual Drop
- E4413** - Lumbar Manual Drop
- E4414** - Pelvic Manual Drop

STANDARD UPHOLSTERY CHOICES:

- | | |
|--|---|
| <input type="checkbox"/>  Imperial Blue | <input type="checkbox"/>  Charcoal |
| <input type="checkbox"/>  Blue Ridge | <input type="checkbox"/>  Mocha |
| <input type="checkbox"/>  Forest | <input type="checkbox"/>  Black |
| <input type="checkbox"/>  Burgundy | <input type="checkbox"/>  Lark |
| <input type="checkbox"/>  Dove | <input type="checkbox"/>  Putty |



E9017