

CORFIT SYSTEM® LS BACK SUPPORT

- Helps relieve lower back pain and discomfort
- Proportional sizing system for a proper fit
- Helps stabilize the spine
- Easy closure assists provide leverage with minimal effort
- Promotes proper posture
- Flexible stays help increase support and stability
- Adjustable elastic pulls for extra compression
- Contours to the waist and hips for comfort



Available From:



HELPS RELIEVE LOWER BACK PAIN

The CorFit System® LS Back Support provides compression to help stabilize the spine and ease lower back pain by creating intra-abdominal pressure. It also helps improve posture and serves as a reminder to use proper lifting techniques by limiting the motion of muscles and ligaments.

PATENTED, PROPORTIONAL SIZING SYSTEM

Our patented, proportional sizing system provides a proper fit in height and length to accommodate various body sizes. Flexible internal stays provide support and stability as well as help prevent the sides from rolling. Easy closure assists provide additional leverage with minimal effort for extra compression and a secure fit; ideal for users with arthritis or weakened muscles. The soft elastic CorEdge® finish on top and bottom edges offers added comfort and flexibility.

COMFORT, COMPRESSION, FLEXIBILITY, & STABILITY

Soft material wicks away moisture, keeping the body cool and comfortable. The anatomical form contours to the natural curves of the waist and hips for comfort. Dual elastic side pulls provide firm compression.

The CorFit system properly fits everyone while providing comfort, compression, flexibility, and stability!

SPECIFICATIONS

Product Selection:
LSB-7000 CorFit System®
LS Back Support

Product Dimensions: XS: 4" x 30",
S: 6" x 35", M/L: 8" x 41", XL: 9" x 46",
2XL: 10" x 52", 3XL: 12" X 58", 4XL: 14" X 64"

Materials: Elastic, heavy loop laminate

Patent: 5,785,672

Latex: This product is not made with natural rubber latex.

Measure around your hips, just below the waistline

SIZE	INCHES	HEIGHT	CM
XS	24-33"	4"	61-84cm
S	27-38"	6"	69-97cm
M/L	32-43"	8"	81-109cm
XL	40-52"	9"	102-132cm
2XL	46-58"	10"	117-147cm
3XL	52-64"	12"	132-163cm
4XL	58-70"	14"	147-178cm