

NAME: _____
COMPANY/ORGANIZATION: _____
DEALER/REP: _____

PHONE: _____
E-MAIL: _____
SALES REP: _____ **QUOTE:** _____



www.phschiropractic.com • 800-743-7738

Thera-P™ Electric Treatment Table

STANDARD FEATURES:


- Load Capacity: 400 lbs.
- 27"W x 77"L x 17-35"H
- Durable and stylish steel base
- White powder coat base
- Whisper quiet lift system
- Height adjustment 17-35"
- Single Foot Control Pedal
- 27" Wide cushions in 6 different configuration options
- 2" High density foam cushions
- Rounded Corners
- Low profile locking casters
- Head Section Angle Range: +40° to -90°
- Lower Extremity Angle Range: 0 to 65°
- Armrest Horizontal Range: 0 to 45°
- Armrest Vertical Range: 0 to 7"
- ADA Compliant

NOTE: International (230V) option available



Ease of use, multiple options and a sleek but durable design make the THERA-P line of electric treatment tables the choice of therapists everywhere. Six top options allow therapists to choose the table that best meets their needs and that ensures patients' most comfortable positioning. The THERA-P electric treatment table is constructed of precision-cut flat and tube steel for durability and style. It also has a whisper-quiet lift system (with elevating foot pedal) and easy-to-operate caster system for better ergonomics.

TABLE CONFIGURATIONS:



- | | | | | | |
|---|---|---|---|---|---|
| <input type="checkbox"/> PT1101
1 Section
27" Width
Qty _____ | <input type="checkbox"/> PT1102
2 Section
27" Width
Qty _____ | <input type="checkbox"/> PT1103
3 Section
27" Width
Qty _____ | <input type="checkbox"/> PT1104
4 Section
27" Width
Qty _____ | <input type="checkbox"/> PT1105
5 Section
27" Width
Qty _____ | <input type="checkbox"/> PT1106
6 Section
27" Width
Qty _____ |
|---|---|---|---|---|---|

OPTIONS:


- PT9027** - 400 lb Lifting Capacity (230V)
- PT9020** - 650 lb Lifting Capacity (115V)
- PT9021** - 650 lb Lifting Capacity (230V)
- PT303x** - 30"W Cushion Configuration
- S1135** - Manual Elevating Midsection
(Only available on PT1103 & PT1105)
- S1077** - Dual Foot Control Pedal
- S1078** - Dual Foot Bar Control



PHSChiropractic
by Pivotal Health Solutions

Thera-P™ Electric Treatment Table

TABLE FEATURES:




Flat Cushion

PT1101



Lower
Extremity
Angle
Ranges from
0 to 65°

PT1102



Lower Extremity
Angle Ranges
from 0 to 65°

Optional
Elevating
Midsection

Head Section
Angle
Ranges from
+40° to -90°


PT1103



Lower Extremity Angle
Ranges from 0 to 65°

Each Head
Section Angle
Ranges from
+40° to -90°

PT1104



Lower Extremity Angle
Ranges from 0 to 65°

Optional Elevating
Midsection

Each Head Section Angle
Ranges from +40° to -90°

Each Armrest Horizontal
Ranges from 0 to 45°

Each Armrest Vertical
Ranges from 0 to 7"

PT1105



Lower Extremity Angle
Ranges from 0 to 65°

Each Head Section Angle
Ranges from +40° to -90°

Each Armrest Vertical
Ranges from 0 to 7"

Each Armrest Horizontal
Ranges from 0 to 45°

PT1106





PHS Chiropractic
by Pivotal Health Solutions

Thera-P™

Electric Treatment Table













STANDARD BASE COLOR

-  White

CUSTOM BASE COLORS:

- | | |
|---|---|
| <input type="checkbox"/>  Black | <input type="checkbox"/>  Green |
| <input type="checkbox"/>  Light Gray | <input type="checkbox"/>  Forest Green |
| <input type="checkbox"/>  Charcoal | <input type="checkbox"/>  Purple |
| <input type="checkbox"/>  Royal Blue | <input type="checkbox"/>  Yellow |
| <input type="checkbox"/>  Navy Blue | <input type="checkbox"/>  Orange |
| <input type="checkbox"/>  Red | <input type="checkbox"/>  Brown |
| <input type="checkbox"/>  Burgundy | |

STANDARD UPHOLSTERY CHOICES:

- | | |
|--|--|
| <input type="checkbox"/>  Royal Blue | <input type="checkbox"/>  American Beauty |
| <input type="checkbox"/>  Imperial Blue | <input type="checkbox"/>  Lark |
| <input type="checkbox"/>  Blue Ridge | <input type="checkbox"/>  Burgundy |
| <input type="checkbox"/>  Amethyst | <input type="checkbox"/>  Dove |
| <input type="checkbox"/>  Forest | <input type="checkbox"/>  Charcoal |
| <input type="checkbox"/>  Black | <input type="checkbox"/>  Mocha |

BRANDING KITS:

- | | | |
|---|---|---|
| <ul style="list-style-type: none"> <input type="checkbox"/> BRANDING KIT 1 - S1170 <ul style="list-style-type: none"> • Slip Cover • Custom Table Base Color • Custom Table Base Graphics <i>*Not available on PT1101, PT1104, or PT1106</i> <input type="checkbox"/> BRANDING KIT 2 - S1171 <ul style="list-style-type: none"> • Slip Cover • White Table Base Color • Custom Table Base Graphics <i>*Not available on PT1101, PT1104, or PT1106</i> | <ul style="list-style-type: none"> <input type="checkbox"/> BRANDING KIT 3 - S1173 <ul style="list-style-type: none"> • Deboss Logo • Custom Table Base Color • Custom Table Base Graphics <i>*Not available on PT1106</i> <input type="checkbox"/> BRANDING KIT 4 - S1174 <ul style="list-style-type: none"> • Deboss Logo • White Table Base Color • Custom Table Base Graphics <i>*Not available on PT1106</i> | <ul style="list-style-type: none"> <input type="checkbox"/> BRANDING KIT 5 - S1176 <ul style="list-style-type: none"> • Color Edge Print • Custom Table Base Color • Custom Table Base Graphics <i>*Not available on PT1106</i> <input type="checkbox"/> BRANDING KIT 6 - S1177 <ul style="list-style-type: none"> • Color Edge Print • White Table Base Color • Custom Table Base Graphics <i>*Not available on PT1106</i> |
|---|---|---|

CUSTOM CUSHION BRANDING:








- PT3027**
Logo Slip Cover
(27"W x 12"H)
- PT3030**
Logo Slip Cover
(30"W x 12"H)
- LOGO DB**
Deboss Logo
(12"W x 12"H)
- DEBPLT12**
Deboss Plate



- LOGO CP**
Color Edge Print
(12"W x 12"H)



COLOR EDGE PRINT UPHOLSTERY CHOICES:

- | | |
|---|--|
| <input type="checkbox"/>  Blue Ribbon | <input type="checkbox"/>  Celestial |
| <input type="checkbox"/>  Purple Passion | <input type="checkbox"/>  Valentine |
| <input type="checkbox"/>  Crimson | <input type="checkbox"/>  Dolphin |
| <input type="checkbox"/>  Black | |

***Color Edge Print upholstery colors do not match Standard Upholstery colors.*

INDIVIDUAL BASE BRANDING:

- S1178**
Table Base Graphics



INDIVIDUAL CUSTOM BASE POWDER COAT:

- PT9026-CC**
Pick Custom Base Color



PT1106 with Dual Footbars